

NITRILE MEDICAL EXAM GLOVES



Nitrile Medical Exam Glove, Powder-Free

- Textured Finish
- Powder-Free
- Latex Free
- Ultra Low Modulus
- Tested for use with Chemotherapy drugs

Certified Model
(NFPA 1999, 2008 Ed)



NitriSkin Blue - MG500 series

SIZES	ORDER NO	PRODUCT SPECIFICATIONS	
XS (5 ½ - 6)	MG500XS	Length min.	9.5" (240mm)
S (6 ½ - 7)	MG500S	Avg. Finger Thickness (mm)	min. 0.15
M (7 ½ - 8)	MG500M	Tensile Strength Before Aging	min. 14 MPa
L (8 ½ - 9)	MG500L	Tensile Strength After Aging	min. 14 MPa
XL (9 ½ - 10)	MG500XL	Elongation Before Aging (%)	min. 500
XXL(10½ - 11)	MG500XXL	Elongation After Aging (%)	min. 500



Characteristics

The NitriSkin Blue glove is a high-performing, latex-free exam glove for use in every application in a healthcare environment. Exceptionally stretchy soft, low modulus nitrile synthetic compound is very comfortable for extended periods. Slightly tacky exterior finish provides excellent gripping properties. Light blue color.

Product Overview

NitriSkin Nitrile Powder Free Exam Gloves offer a high valued and dependable solution to those allergic to natural latex protein. Meets and exceeds the standard emergency glove medical of NFPA 1999 Standards on protective clothing for Emergency medical operations, 2008 Edition. Tested approved by SEI for water leak test that meets all the requirements of ASTM 6319-00 and NFPA 1999, Edition 2008. Tested for use with Chemotherapy drugs.

Features and Benefits

Much softer and stretchy than other nitrile gloves without any latex additives. 100% latex free eliminates potential allergic reactions to natural latex protein. Textured surface provides strong grip for both wet and dry applications. Superior puncture, tear and chemical resistance with excellent tactile-sensitivity and dexterity.

Quality Standards:

- Meets or exceeds the following standards: ASTM D6319-00 (USA), EN 455 (EEC), AS 40 (Australia), FDA, NFPA, UL
- Quality sampled in accordance with MIL.STD 105E
- ISO 9002 Certified Manufacturing

Packaging: 100 gloves per box, 10 boxes per case. (XL & XXL - 95 gloves per box)

Sizes: X-Small, Small, Medium, Large, Extra Large and 2 Extra Large

Storage

Store in original packing in a cool, dry and well-ventilated area, away from dust, sunlight, fluorescent lighting, moisture, X-ray, and excessive heat above 100°F (37°C).

Meets or exceeds ASTM D6319