


NITRILE MEDICAL EXAM GLOVES

NitriSkin® XP Exam Gloves are made from 100% soft pliable nitrile, a synthetic compound that's completely free of skin-irritating rubber proteins. These stretchy gloves conform to your hands far more comfortably than other non-latex alternatives, such as vinyl. So you get a cool, dry, precise fit, with less hand fatigue. NitriSkin® XP Exam Gloves are also more resistant to punctures and tears than either latex or vinyl and stand up to a variety of splashed chemicals, solvents, oils, and bases.

FEATURES

- 8 mil thickness
- 12" Extended Cuff
- Thin Wall nitrile technology
- More puncture-resistant
- Superior Strength and durability
- Textured finish for superior grip
- Soft and Stretchy
- Excellent alternative to NBR
- Higher resistance to certain chemicals or solvents
- Dated Lot Codes for quality assurance and traceability
- Safety Sealed Packaging
- Tested for use with Chemotherapy drugs.



Certified Model
(NFPA 1999, 2008 Ed)



NitriSkin XP - MG5008 series

SIZES

- S (6 ½ - 7)
- M (7 ½ - 8)
- L (8 ½ - 9)
- XL (9 ½ - 10)
- 2XL (10 ½ - 11)
- 3XL (11 ½ - 12)

ORDER NO

- MG5008S
- MG5008M
- MG5008L
- MG5008XL
- MG5008XXL
- MG5008XXXL

PRODUCT SPECIFICATIONS

Length min.	12"
Finger Thickness (mm/mil)	0.21 / 8.2
Cuff Thickness (mm/mil)	0.127 / 5
Palm Thickness (mm/mil)	0.152 / 6
Tensile Strength Before Aging	min. 14 MPa
Tensile Strength After Aging	min. 14 MPa
Elongation Before Aging (%)	min. 500
Elongation After Aging (%)	min. 500

Characteristics

When a high risk situation requires the highest standard in protection. These nitrile synthetic gloves are typically used for tasks involving prolonged exposure to blood, body fluids, chemotherapeutic agents, surgical instrument cleaning solutions and other chemicals. Excellent donning properties and a textured, slightly tacky surface for a better grip in wet or dry conditions.

Ingredients:

100% NBR Nitrile butadiene synthetic rubber formulation

Product Overview

NitriSkin Nitrile Powder Free Exam Gloves offer a high valued and dependable solution to those allergic to natural latex protein. Meets and exceeds the standard emergency glove medical of NFPA 1999 Standards on protective clothing for Emergency medical operations, 2008 Edition. Tested approved by by SEI for water leak test that meets all the requirements of ASTM 6319-00 and NFPA 1999, Edition 2008. Tested for use with Chemotherapy drugs.

Quality Standards:

- Meets or exceeds the following standards: ASTM D6319-00 (USA), EN 455 (EEC), AS 40 (Australia), FDA, NFPA, UL
- Quality sampled in accordance with MIL.STD 105E
- ISO 9002 Certified Manufacturing

Packaging:

50 gloves per box, 10 boxes per case.

Sizes:

Small, Medium, Large, Extra Large and Extra Extra Large