


**STERILE LATEX-FREE EXAM GLOVES**


### NitriSkin\* Plus Powder-Free Nitrile Surgical Glove

- Textured Finish
- Tactile Sensivity compared to Latex
- Latex-Free
- Excellent Barrier

Cerified Model  
(NFPA 1999, 2008 Ed)



#### Packaging

25 pairs per box, 4 boxes per case.

#### Sizes

6.0 to 9.0

### NitriSkin Plus PF Surgical - MGS50 series

SIZES	ORDER #	PRODUCT SPECIFICATIONS	
6.0	MGS5060	Length	min. 11½" (290mm)
6.5	MGS5065	Avg. Finger Thickness (mm)	min. 0.22
7.0	MGS5070	Tensile Strength Before Aging	min. 24 MPa
7.5	MGS5075	Tensile Strength After Aging	min. 18 MPa
8.0	MGS5080	Elongation Before Aging (%)	min. 800
8.5	MGS5085	Elongation After Aging (%)	min. 760
9.0	MGS5090		

#### Characteristics

The NitriSkin Plus Nitrile Surgical gloves contain no natural rubber latex proteins. It's soft, highly elastic formulation provides the maximum in comfort, fit, and protection, while offering exceptional tactile sensitivity. Hand Specific. Textured Finish.

#### Features and Benefits

Nitrile Surgical Gloves are intended for use in environments within hospitals and other healthcare facilities. The gloves are appropriate for use during invasive and non-invasive medical procedures requiring sterility. They are intended to be worn by operating room personnel to protect a surgical wound from contamination.

#### Storage

Store in original packing in a cool, dry and well-ventilated area, away from dust, sunlight, fluorescent lighting, moisture, X-ray, and excessive heat above 100°F (37°C).