



STERILE EXAM GLOVES



Sterile Singles, Powder-Free Nitrile Exam

SIZES	REORDER NO.	PRODUCT SPECIFICATIONS
S (6 ½ - 7)	MGSE502S	Length min. 9.5" (240mm)
M (7 ½ - 8)	MGSE502M	Avg. Finger Thickness (mm) min. 0.20
L (8 ½ - 9)	MGSE502L	Tensile Strength Before Aging min. 28 MPa
XL (9 ½ - 10)	MGSE502XL	Tensile Strength After Aging min. 20 MPa
		Elongation Before Aging (%) min. 900
		Elongation After Aging (%) min. 880

NitriSkin Sterile Singles

- Sterile, Powder-Free
- Ambidextrous
- Latex-free

**Packaging :** 100 gloves per box, 4 boxes per case.  
**Sizes:** Small, Medium, Large and Extra Large



Sterile Pairs, Powder-Free Nitrile Exam

SIZES	REORDER NO.	PRODUCT SPECIFICATIONS
S (6 ½ - 7)	MGSE552S	Length min. 9.5" (240mm)
M (7 ½ - 8)	MGSE552M	Avg. Finger Thickness (mm) min. 0.20
L (8 ½ - 9)	MGSE552L	Tensile Strength Before Aging min. 28 MPa
XL (9 ½ - 10)	MGSE552XL	Tensile Strength After Aging min. 20 MPa
		Elongation Before Aging (%) min. 900
		Elongation After Aging (%) min. 880

NitriSkin Sterile Pairs

- Sterile, Powder-Free
- Ambidextrous
- Latex-free

**Packaging :** 50 pairs per box, 4 boxes per case.  
**Sizes:** Small, Medium, Large and Extra Large



Product Overview

NitriSkin Sterile Singles and Pairs are now the ultimate choice use of medical professionals in high demand environments like Obstetricians and Gynecologists, intensive care unit ICC/CCU, traumas, burn centers and wound care that requires a high quality latex-free sterile glove that delivers maximum barrier protection.

Features and Benefits

Powder Free minimizes the risk of developing chemical residue related allergies. Eliminate powder related contamination in critical exam or in clean room environments. Our ambidextrous sterile latex-free exam gloves features excellent tactile sensitivity with finger textured to provide a sure grip.

Tested for use with Chemotherapy drugs\*

\*\* Additional information available upon requests