


**NITRILE MEDICAL EXAM GLOVES**

**Neutron Nitrile Powder-Free Exam Glove**

- Textured Finish
- Tactile Sensivity compared to Latex
- Latex-Free
- Excellent Barrier


**Neutron Nitrile - MNE505 series**

| SIZES         | ORDER NO | PRODUCT SPECIFICATIONS        |              |
|---------------|----------|-------------------------------|--------------|
| XS (5 ½ - 6)  | MNE5051  | Length min.                   | 9.5" (240mm) |
| S (6 ½ - 7)   | MNE5052  | Avg. Finger Thickness (mm)    | min. 0.09    |
| M (7 ½ - 8)   | MNE5053  | Tensile Strength Before Aging | min. 18 MPa  |
| L (8 ½ - 9)   | MNE5054  | Tensile Strength After Aging  | min. 16 MPa  |
| XL (9 ½ - 10) | MNE5055  | Elongation Before Aging (%)   | min. 600     |
| XXL(10½ - 11) | MNE5056  | Elongation After Aging (%)    | min. 550     |

**Characteristics**

Medgluv new latex-free Exam, Neutron Nitrile will offer you the protection of Nitrile and feeling of latex sensitivity. Freatured with an advanced stretch fomulation, it will provide you ultimate comfort for extended periods. Packed in 200 gloves per box, it will increase productivity and maximize available storage space.

**Product Overview**

This new Latex-Free Exam Gloves offer a high valued and dependable solution to those allergic to natural latex protein. These gloves, formulated with 100% synthetic nitrile polymer, are latex free, powder free and odor free. The unique thermal sensitive polymer compound absorbs body heat and easily conforms to user's hands; thus, exerts virtually no pressure on the hands and reduces fatigue after extended wear.

**Features and Benefits**

Much softer and elastic than other nitrile gloves without any latex additives. 100% latex free eliminates potential allergic reactions to natural latex protein. Textured surface provides strong grip for both wet and dry applications. Superior puncture, tear and chemical resistance with excellent tactile-sensitivity and dexterity.

**Packaging**

100 gloves per box, 10 boxes per case. (XL & XXL - 100 gloves per box)

**Sizes**

X-Small, Small, Medium, Large, Extra Large and 2 Extra Large

**Storage**

Store in original packing in a cool, dry and well-ventilated area, away from dust, sunlight, fluorescent lighting, moisture, X-ray, and excessive heat above 100°F (37°C).

**THIS GLOVE MEETS THE EMERGENCY MEDICAL GLOVE REQUIREMENTS OF NFPA 1999, STANDARD ON PROTECTIVE CLOTHING FOR EMERGENCY MEDICAL OPERATIONS, 2008 EDITION. THIS GLOVE IS FOR SINGLE USE ONLY**