



## LATEX MEDICAL EXAM GLOVES



Plus Latex Exam, Powder-Free

- Textured Finish
- Powder-Free
- Traditional Chlorination
- Low Protein (<50µg\*)



### NeuGrip PF - MG105 Series

SIZES	REORDER NO.	PRODUCT SPECIFICATIONS	
XS (5 ½ - 6)	MG105XS	Length	min. 9.5" (240mm)
S (6 ½ - 7)	MG105S	Avg. Finger Thickness (mm)	min. 0.19
M (7 ½ - 8)	MG105M	Tensile Strength Before Aging	min. 24 MPa
L (8 ½ - 9)	MG105L	Tensile Strength After Aging	min. 18 MPa
XL (9 ½ - 10)	MG105XL	Elongation Before Aging (%)	min. 850
		Elongation After Aging (%)	min. 820

#### Characteristics

NeuGrip powder-free product. Fully textured to enhance grip. Durable, medium modulus latex is 100% inspected for consistency. A high-quality, yet economical choice

#### Product Overview

NeuGrip Latex Powder Free Exam Gloves are the supreme choice of medical professionals who demand uncompromising quality and protection, with the minimum risk of potential glove related chemical or protein allergies. Specially manufactured with light chlorination rinse, these gloves maintain a non-slip, 'tacky' surface that provides excellent grip for both dry and wet applications.

Our NeuGrip Powder-Free Latex Exam glove is an early recipient of the prestigious **SMG** certification from the Malaysian Rubber Research Institute as recognition for its compliance with high quality standards.

#### Features and Benefits

Powder Free minimizes the risk of developing latex protein or chemical residue related allergies. Eliminate powder related contamination in critical exam or in clean room environment. Full textured surface provides superior grip in wet or dry applications; Smooth interior finish for easy donning and comfortable wear.

**Packaging:** 100 gloves per box, 10 boxes per case.

**Sizes:** X-Small, Small, Medium, Large and Extra Large

**Caution:** *No sensitizing threshold to natural rubber latex has been established. Any individual with hypersensitivity to natural rubber latex should consult a physician before using any latex product.*