

LATEX MEDICAL EXAM GLOVES



TuffSkin XP Latex Exam, Extended Cuff

- Textured Finish
- Powder-Free
- 12", 15 mil
- Low Protein (<50µg*)



TuffSkin XP - MG1215 Series

SIZES	REORDER NO.	PRODUCT SPECIFICATIONS	
S (6½ - 7)	MG1215S	Length	min. 12" (305 mm)
M (7½ - 8)	MG1215M	Avg. Finger Thickness (mm)	min. 0.36
L (8½ - 9)	MG1215L	Tensile Strength Before Aging	min. 24 MPa
XL (9½ - 10)	MG1215XL	Tensile Strength After Aging	min. 20 MPa
2XL (10½ - 11)	MG1215XXL	Elongation Before Aging (%)	min. 750
		Elongation After Aging (%)	min. 600

Characteristics

TuffSkin XP Powder-Free Latex Exam glove is 4 times thicker than standard exam gloves, and offers added durability in the palm and fingers. An extended cuff offers greater coverage over the wrist and arm. TuffSkin XP Glove meets or exceeds the emergency medical glove requirements of NFPA 1999, 2008 edition.

Product Overview

For maximum protection in high risk environments. Chlorinated and coated for rapid donning in critical situations. Textured finish with a medium surface tack for positive grip. High tensile strength for added protection. Our Premium Powder-Free Latex Exam glove is an early recipient of the prestigious **SMG** certification from the Malaysian Rubber Research Institute as recognition for its compliance with high quality standards.

Features and Benefits

Powder Free minimizes the risk of developing latex protein or chemical residue related allergies. Eliminate powder related contamination in critical exam or in clean room environment. Full textured surface provides superior grip in wet or dry applications. Extended cuff and 15 mil thickness for maximum high risk protection.

Packaging

50 gloves per box, 10 boxes per case.

Sizes

X-Small, Small, Medium, Large, Extra Large, 2 Extra Large

Quality Standards:

- Meets or exceeds the following standards: ASTM D35780 (USA), EN 455 (EEC), AS 40 (Australia), FDA, NFPA, UL
- Quality sampled in accordance with MIL.STD 105E
- ISO 9002 Certified Manufacturing



THIS GLOVE MEETS THE EMERGENCY MEDICAL GLOVE REQUIREMENTS OF NFPA 1999, STANDARD ON PROTECTIVE CLOTHING FOR EMERGENCY MEDICAL OPERATIONS, 2008 EDITION. THIS GLOVE IS FOR SINGLE USE ONLY

CAUTION: This product contains natural rubber latex which may cause allergic reactions. Safe use of these gloves by or on latex sensitized individuals has not been established.