

SYNTHETIC VINYL MEDICAL EXAM GLOVES



NeuSkin Vinyl Lightly Powdered

- Smooth Finish
- Lightly Powdered
- Latex Free

NeuSkin Vinyl Powdered - MG350 series

SIZES	ORDER NO	PRODUCT SPECIFICATIONS	
XS (5 ½ - 6)	MG350XS	Length	min. 9.5" (240mm)
S (6 ½ - 7)	MG350S	Avg. Finger Thickness (mm)	min. 0.15
M (7 ½ - 8)	MG350M	Tensile Strength Before Aging	min. 10 MPa
L (8 ½ - 9)	MG350L	Tensile Strength After Aging	min. 10 MPa
XL (9 ½ - 10)	MG350XL	Elongation Before Aging (%)	min. 400
		Elongation After Aging (%)	min. 400

Characteristics

Synthetic material contains no natural rubber latex proteins. Pure PVC polymer is soft and semi-elastic. Light Powder for easy donning. Recommended low cost choice for non-latex requirements.

Product Overview

Developed from a special PVC material to feature exceptional softness and dexterity for a relaxed and comfortable fit, NeuSkin Vinyl Lightly Powdered Exam gloves offer an affordable alternative to individuals who prefer non-latex protection.

Features and Benefits

Powdered with USP grade edible cornstarch for easy donning.
100% latex free eliminates the risk of developing latex protein related allergies.
Ergonomically designed to fit well for maximum comfort. Very soft yet also very durable.

Packaging

100 gloves per box, 10 boxes per case.

Sizes

X-Small, Small, Medium, Large and Extra Large

Storage

Store in original packing in a cool, dry and well-ventilated area, away from dust, sunlight, fluorescent lighting, moisture, X-ray, and excessive heat above 100°F (37°C).