


STERILE EXAM GLOVES



Sterile Singles, Powder-Free Latex Exam

- Sterile, Powder-Free
- Ambidextrous
- Low Protein (<50µg*)

Packaging

100 gloves per box, 4 boxes per case.

Sizes

Small, Medium, Large and Extra Large

Sterigluv Sterile Singles, PF Latex

SIZES

- S (6 ½ - 7)
- M (7 ½ - 8)
- L (8 ½ - 9)
- XL (9 ½ - 10)

REORDER NO.

- MGSE100S
- MGSE100M
- MGSE100L
- MGSE100XL

PRODUCT SPECIFICATIONS

| | |
|-------------------------------|-------------------|
| Length | min. 9.5" (240mm) |
| Avg. Finger Thickness (mm) | min. 0.20 |
| Tensile Strength Before Aging | min. 28 MPa |
| Tensile Strength After Aging | min. 20 MPa |
| Elongation Before Aging (%) | min. 900 |
| Elongation After Aging (%) | min. 880 |

Product Overview

Sterigluv* Sterile Singles are now the ultimate choice use of medical professionals in high demand environments like Obstetricians and Gynecologists, intensive care unit ICC/CCU, traumas, burn centers and wound care that requires a high quality sterile glove that delivers maximum barrier protection.

Features and Benefits

Powder Free minimizes the risk of developing latex protein or chemical residue related allergies. Eliminate powder related contamination in critical exam or in clean room environment. Our ambidextrous sterile latex exam gloves feature excellent tactile sensitivity to provide a sure grip.

Caution: *No sensitizing threshold to natural rubber latex has been established. Any individual with hypersensitivity to natural rubber latex should consult a physician before using any latex product.*