

NitriSkin XP



Extended Cuff PF Nitrile Exam Gloves

NitriSkin XP Nitrile Exam Gloves is a thicker quality nitrile glove that resists punctures and tears to provide you outstanding barrier protection that is full textured for a more confident grip. When your daily job exposes you to blood-borne pathogens, bacteria and fluids, NitriSkin XP is a glove you can trust.

MAXIMUM BARRIER PROTECTION FOR ALL RISKS

Description:

- When a high-risk situation requires the highest standard in protection.
- Offer excellent sensitivity for touch perception and dependable protection
- Meet standards of the National Fire Protection Association (NFPA)
- Meets and exceeds NIOSH recommendations for Fentanyl.
- Meets both ASTM D6978-05 and USP 800 Personal Protection Equipment requirement for use of chemotherapy gloves when handling hazardous drugs.



Item	Description	Size	Packaging
MG5008S	NitriSkin XP Nitrile Exam Glove	Small	50 eaches/box 10 boxes/case
MG5008M	NitriSkin XP Nitrile Exam Glove	Medium	50 eaches/box 10 boxes/case
MG5008L	NitriSkin XP Nitrile Exam Glove	Large	50 eaches/box 10 boxes/case
MG5008XL	NitriSkin XP Nitrile Exam Glove	X-Large	50 eaches/box 10 boxes/case
MG5008XXL	NitriSkin XP Nitrile Exam Glove	2X-Large	45 eaches/box 10 boxes/case
MG5008XXXL	NitriSkin XP Nitrile Exam Glove	3X-Large	45 eaches/box 10 boxes/case

Product Specifications

Gauge Thickness	MM	MIL
Middle Finger:	.21	8.2
Palm:	.17	6.8
Cuff:	.14	5.5
Average Length	305mm	12"

Physical Properties

Before Aging	
Tensile Strength:	35 MPa
Ultimate Elongation:	600%
After Aging	
Tensile Strength:	32 MPa
Ultimate Elongation:	550%

Quality Standards

Meet and exceeds ASTM D6319 and F1671: Standard Specification for Nitrile: Examination Gloves for Medical Application and Standard Test method for Resistance of Materials used in Protective Clothing in Penetration by Blood-Borne Pathogens Using Phi-X174 Bacteriophage as a Test System.

Note: Specifications are subject to change without notice.



Tested with 30 CHEMO DRUGS
in accordance with ASTM D6978 and barrier tested against Fentanyl exposure.



This glove meets **NFPA 1999**, Standard on Protective Clothing and Ensembles for Emergency Medical Operations.

For additional information, please visit: www.medgluv.com
Call 1-866-MEDGLUV (1-866-633-4588).

Keep out of sunlight. Store in a cool dry place. Keep away from sources of ozone and ignition.



www.medgluv.com

