

NitriSkin XP



Extended Cuff PF Nitrile Exam Gloves

NitriSkin XP Nitrile Exam Gloves is a thicker quality nitrile glove that resists punctures and tears to provide you outstanding barrier protection that is full textured for a more confident grip. When your daily job exposes you to blood-borne pathogens, bacteria and fluids, NitriSkin XP is a glove you can trust

Description:

- When a high-risk situation requires the highest standard in protection.
- Offer excellent sensitivity for touch perception and dependable protection
- Meet standards of the National Fire Protection Association (NFPA)
- Meets and exceeds NIOSH recommendations for Fentanyl.
- Meets both ASTM D6978-05 and USP 800 Personal Protection Equipment requirement for use of chemotherapy gloves when handling hazardous drugs.

Item	Description	Size	Packaging
MG5008XS	NitriSkin XP Nitrile Exam Glove	X-Small	50 eaches/box 10 boxes/case
MG5008S	NitriSkin XP Nitrile Exam Glove	Small	50 eaches/box 10 boxes/case
MG5008M	NitriSkin XP Nitrile Exam Glove	Medium	50 eaches/box 10 boxes/case
MG5008L	NitriSkin XP Nitrile Exam Glove	Large	50 eaches/box 10 boxes/case
MG5008XL	NitriSkin XP Nitrile Exam Glove	X-Large	50 eaches/box 10 boxes/case
MG5008XXL	NitriSkin XP Nitrile Exam Glove	2X-Large	45 eaches/box 10 boxes/case
MG5008XXXL	NitriSkin XP Nitrile Exam Glove	3X-Large	45 eaches/box 10 boxes/case

Product Specifications

Gauge Thickness	MM	MIL
Middle Finger:	.21	8.2
Palm:	.17	6.8
Cuff:	.14	5.5
Average Length	305mm	12"

Physical Properties

Before Aging	
Tensile Strength:	35 MPa
Ultimate Elongation:	600%
After Aging	
Tensile Strength:	32 MPa
Ultimate Elongation:	550%

Quality Standards

Meet and exceeds ASTM D6319 and F1671: Standard Specification for Nitrile: Examination Gloves for Medical Application and Standard Test method for Resistance of Materials used in Protective Clothing in Penetration by Blood-Borne Pathogens Using Phi-X174 Bacteriophage as a Test System.

Note: Specifications are subject to change without notice.



www.medgluv.com

MAXIMUM BARRIER PROTECTION FOR ALL RISKS



Tested with 53 CHEMO DRUGS
in accordance with ASTM D6978 and barrier tested against Fentanyl exposure.



This glove meets **NFPA 1999**, Standard on Protective Clothing and Ensembles for Emergency Medical Operations.

For additional information, please visit: www.medgluv.com
Call 1-866-MEDGLUV (1-866-633-4588).

Keep out of sunlight. Store in a cool dry place. Keep away from sources of ozone and ignition.

