

## 2-Ply Blue/White pH Nitrile Exam Gloves

These 2-Ply technology outer Blue and inner White elongated latex-free nitrile gloves offers the unique ability to detect tears and rips by showing the high visibility white inner through the breach.

### Description:

- Designed specially for the Emergency Medical Profession for high risk situation that requires the highest standards in double protection and safety in administering critical and patient care.
- Breach alert dual-layer provides extra durability and safety.
- Ideal for EMS, Law Enforcement and Corrections
- Meets and exceeds NIOSH recommendations for Fentanyl and gastric acid resistant.
- Meets both ASTM D6978-05 and USP 800 Personal Protection Equipment requirement for use of chemotherapy gloves when handling hazardous drugs.
- Meets and exceeds the standard emergency glove medical of NFPA 1999 Standards on protective clothing, Emergency Medical Operations, Edition 2018

| Item    | Description                             | Size     | Packaging                      |
|---------|---|----------|--------------------------------|
| MG50091 | NitraGrip Pro, 2-Ply Nitrile Exam Glove | Small    | 50 eaches/box<br>10 boxes/case |
| MG50092 | NitraGrip Pro, 2-Ply Nitrile Exam Glove | Medium   | 50 eaches/box<br>10 boxes/case |
| MG50093 | NitraGrip Pro, 2-Ply Nitrile Exam Glove | Large    | 50 eaches/box<br>10 boxes/case |
| MG50094 | NitraGrip Pro, 2-Ply Nitrile Exam Glove | X-Large  | 50 eaches/box<br>10 boxes/case |
| MG50095 | NitraGrip Pro, 2-Ply Nitrile Exam Glove | 2X-Large | 45 eaches/box<br>10 boxes/case |
| MG50096 | NitraGrip Pro, 2-Ply Nitrile Exam Glove | 3X-Large | 45eaches/box<br>10 boxes/case  |

### Product Specifications

| Gauge Thickness | MM    | MIL |
|-----------------|-------|-----|
| Middle Finger:  | .21   | 8.2 |
| Palm:           | .17   | 6.8 |
| Cuff:           | .14   | 5.5 |
| Average Length  | 305mm | 12" |

### Physical Properties

| Before Aging         |        |
|----------------------|--------|
| Tensile Strength:    | 35 MPa |
| Ultimate Elongation: | 600%   |
| After Aging          |        |
| Tensile Strength:    | 32 MPa |
| Ultimate Elongation: | 550%   |

### Quality Standards

Meet and exceeds ASTM D6319 and F1671: Standard Specification for Nitrile: Examination Gloves for Medical Application and Standard Test method for Resistance of Materials used in Protective Clothing in Penetration by Blood-Borne Pathogens Using Phi-X174 Bacteriophage as a Test System.

Note: Specifications are subject to change without notice.



## EXCELLENT BARRIER PROTECTION FOR ALL RISKS



**Tested with 53 CHEMO DRUGS**  
in accordance with ASTM D6978 and barrier tested against Fentanyl exposure.



This glove meets **NFPA 1999**, Standard on Protective Clothing and Ensembles for Emergency Medical Operations.

For additional information, please visit: [www.medgluv.com](http://www.medgluv.com)  
Call 1-866-MEDGLUV (1-866-633-4588).

Keep out of sunlight. Store in a cool dry place. Keep away from sources of ozone and ignition.

