

# NitraPro 16"



## Elongated Nitrile Exam Gloves

Excellent for sterilization, decontamination, disinfect, cleaning or handling with chemicals in a hazardous environment. This 16-inch specialty glove is designed to offer exceptional added coverage and protection where there is potential exposure to sharps and chemicals. Tested with a wide range of chemotherapy drugs including Carmustine and Thiotepea.

### Description:

- Full textured fingertips enhance grip in wet and dry situations.
- Longer 16" extended length
- Tested for puncture resistance and color-coded sizing on box.
- Meets and exceeds NIOSH recommendations for Fentanyl and gastric acid resistant.
- Meets both ASTM D6978-05 and USP 800 Personal Protection Equipment requirement for use of chemotherapy gloves when handling hazardous drugs.

Item	Description	Size	Packaging
MG50161	NitraPro 16 Nitrile Exam Glove	Small	50 eaches/box 10 boxes/case
MG50162	NitraPro 16 Nitrile Exam Glove	Medium	50 eaches/box 10 boxes/case
MG50163	NitraPro 16 Nitrile Exam Glove	Large	50 eaches/box 10 boxes/case
MG50164	NitraPro 16 Nitrile Exam Glove	X-Large	50 eaches/box 10 boxes/case

### Product Specifications

Gauge Thickness	MM	MIL
Middle Finger:	0.29	11.5
Palm:	0.255	10.1
Cuff:	0.22	8.6
Average Length	405mm	16"

### Physical Properties

Before Aging	
Tensile Strength:	32 MPa
Ultimate Elongation:	650%
After Aging	
Tensile Strength:	25 MPa
Ultimate Elongation:	600%

### Quality Standards

Manufactured in compliance with multiple international standards, including ASTM D6319, ASTM F1671, ASTM D5151 and ASTM D6978

Note: Specifications are subject to change without notice.

## MAXIMUM BARRIER PROTECTION FOR ALL RISKS



view this product  
on our website:



**Tested with 24 CHEMO DRUGS**  
in accordance with ASTM D6978 and barrier tested against Fentanyl exposure.

For additional information, please visit: [www.medgluv.com](http://www.medgluv.com)  
Call 1-866-MEDGLUV (1-866-633-4588).

Keep out of sunlight. Store in a cool dry place. Keep away from sources of ozone and ignition.



[www.medgluv.com](http://www.medgluv.com)

