

# Dragon Skinz

## Powder-Free Nitrile Exam Gloves



Dragon Skinz nitrile glove is a high performing, latex-free exam glove for use in every application in every environment. Exceptionally stretchy soft, low modulus nitrile synthetic compound excellent for extended periods. Slightly tacky exterior finish provides excellent gripping properties.

### Description:

- Dragon Skinz black exam gloves are strong, highly visible gloves provide the advanced protection and optimal grip needed for demanding, high-stress situations.
- Powder-Free nitrile gloves offer dependable strength and are fully textured for an enhanced grip in demanding environments.
- Meets and exceeds NIOSH recommendations for Fentanyl.
- Meets both ASTM D6978-05 and USP 800 Personal Protection Equipment requirement for use of chemotherapy gloves when handling hazardous drugs.

| Item   | Description                     | Size     | Packaging                       |
|--------|---------------------------------|----------|---------------------------------|
| MG5990 | Dragon Skinz Nitrile Exam Glove | X-Small  | 100 eaches/box<br>10 boxes/case |
| MG5991 | Dragon Skinz Nitrile Exam Glove | Small    | 100 eaches/box<br>10 boxes/case |
| MG5992 | Dragon Skinz Nitrile Exam Glove | Medium   | 100 eaches/box<br>10 boxes/case |
| MG5993 | Dragon Skinz Nitrile Exam Glove | Large    | 100 eaches/box<br>10 boxes/case |
| MG5994 | Dragon Skinz Nitrile Exam Glove | X-Large  | 95 eaches/box<br>10 boxes/case  |
| MG5995 | Dragon Skinz Nitrile Exam Glove | 2X-Large | 90 eaches/box<br>10 boxes/case  |

## MAXIMUM BARRIER PROTECTION FOR ALL RISKS



### Product Specifications

| Gauge Thickness | MM    | MIL  |
|-----------------|-------|------|
| Middle Finger:  | .14   | 5.5  |
| Palm:           | .12   | 4.7  |
| Cuff:           | .09   | 3.5  |
| Average Length  | 242mm | 9.5" |

### Physical Properties

#### Before Aging

|                      |        |
|----------------------|--------|
| Tensile Strength:    | 38 MPa |
| Ultimate Elongation: | 550%   |

#### After Aging

|                      |        |
|----------------------|--------|
| Tensile Strength:    | 33 MPa |
| Ultimate Elongation: | 500%   |

### Quality Standards

THIS GLOVE IS FOR SINGLE USE ONLY. THIS GLOVE MEETS THE SINGLE-USE EMERGENCY MEDICAL EXAMINATION GLOVE REQUIREMENTS OF NFPA 1999, STANDARD ON PROTECTIVE CLOTHING AND ENSEMBLES FOR EMERGENCY MEDICAL OPERATIONS, 2018 EDITION.

Note: Specifications are subject to change without notice.



National Fire  
Protection Agency.



**Tested with 30 CHEMO DRUGS**  
in accordance with ASTM D6978 and barrier  
tested against Fentanyl exposure.

For additional information, please visit: [www.medgluv.com](http://www.medgluv.com)  
Call 1-866-MEDGLUV (1-866-633-4588).

Keep out of sunlight. Store in a cool dry place. Keep away from sources of ozone and ignition.



[www.medgluv.com](http://www.medgluv.com)

