

NitraGrip



Nitrile Sterile Exam Pairs Extended Cuff

NitraGrip Sterile Gloves is made from synthetic nitrile formulation that provides high elasticity with superior barrier protection and puncture resistance. It is an ideal solution for a range of sterile examination (non-surgical) procedures.

Description:

- Perfect for USP 797 and 800 Standards Compliance for pharmacy compounding and sterile processing.
- Sterile procedures outside the OR
- Meets and exceed NIOSH recommendations for Fentanyl.
- Meets both ASTM D6978-05 and USP 800 Personal Protection Equipment requirement for use of chemotherapy gloves when handling hazardous drugs.

Item	Description	Size	Packaging
MGSE5121	NitraGrip Nitrile Sterile Pairs Ext Cuff	Small	50 pairs/box 4 boxes/case
MGSE5122	NitraGrip Nitrile Sterile Pairs Ext Cuff	Medium	50 pairs/box 4 boxes/case
MGSE5123	NitraGrip Nitrile Sterile Pairs Ext Cuff	Large	50 pairs/box 4 boxes/case
MGSE5124	NitraGrip Nitrile Sterile Pairs Ext Cuff	X-Large	50 pairs/box 4 boxes/case



Product Specifications

Gauge Thickness	MM
Middle Finger:	6.3 mil
Palm:	4.5 mil
Cuff:	.10
Average Length	305mm

Physical Properties


Before Aging	
Tensile Strength:	35 MPa
Ultimate Elongation:	600%
After Aging	
Tensile Strength:	33 MPa
Ultimate Elongation:	560%

Quality Standards

Meet and exceeds ASTM D6319 and F1671: Standard Specification for Nitrile: Examination Gloves for Medical Application and Standard Test method for Resistance of Materials used in Protective Clothing in Penetration by Blood-Borne Pathogens Using Phi-X174 Bacteriophage as a Test System.

Note: Specifications are subject to change without notice.



 **Tested with 53 CHEMO DRUGS**
in accordance with ASTM D6978 and barrier tested against Fentanyl exposure.

For additional information, please visit: www.medgluv.com
Call 1-866-MEDGLUV (1-866-633-4588).

Keep out of sunlight. Store in a cool dry place. Keep away from sources of ozone and ignition.



www.medgluv.com

