

## AAMI Level 1 Isolation Gown - MGGL104

Disposable Protective Isolation Gown. Non-Sterile 45 grams Polyethylene film fluid-resistance protective gown. Impervious to the penetration of Fluid liquids.

- Disposable non-surgical isolation gown
- Breathable, Knitted cuff wrists
- Straps on neck and waists
- Impervious to Fluids.
- 100 gowns per case.



## AAMI Level 2 Isolation Gown - MGGL24


Disposable Protective Isolation Gown, non-sterile, 27 grams polypropylene non-woven fabric. Impervious to the penetration of fluid liquids.

- Disposable non-surgical isolation gown
- Breathable, Knitted cuff wrists
- Straps on neck and waists
- Impervious to Fluids.
- Ultrasonic sewn seams
- Meets and Exceeds AAMI PB70 – Level 2
- 100 gowns per case.

## AAMI Level 3 Isolation Gown - MGGL34

Disposable Protective Isolation Gown, non-sterile, 45 grams polypropylene with polyethylene coating. Ultrasonic sewn seams. Impervious to the penetration of fluid liquids.

- Disposable Non-Surgical Isolation Gown
- 100 gowns per case
- Splash Resistance
- Straps on neck and waist
- Elastic Knit Rib Cuff
- Ultrasonic sewn seams.



Certified Gowns Test Results:

AAMI Level	Test	Liquid Challenge	Result	Medgluv Gown Test Result
1	AATCC 42 Impact Penetration	Water	≤ 4.5 g	0.06 g PASS
2	AATCC 42 Impact Penetration	Water	≤ 1.0 g	0.06 g PASS
	AATCC 127 Hydrostatic Pressure	Water	≤ 20 cm	144.97 cm PASS
3	AATCC 42 Impact Penetration	Water	≤ 1.0 g	0.06 g PASS
	AATCC 127 Hydrostatic Pressure	Water	< 5.0 cm	144.49 cm PASS

