



Powder-Free Latex Exam Gloves

NeuGrip Latex Aloe Powder-Free Exam Gloves have excellent barrier properties combined with great fit, feel and are powder-free to minimize skin irritation. They offer a soft and comfortable feel with outstanding tactile sensitivity. Naturally helps moisturize and soothe dry, chapped hands.

**EXCEPTIONAL FIT,
THICKNESS FOR
DURABLE PROTECTION**

Description:

- Top quality Latex Exam glove with the benefits of aloe vera
- Naturally helps moisturize and soothe dry, chapped hands
- Strong durable powder-free natural rubber latex gloves provide superior barrier protection to minimize skin irritation.
- Standard Malaysian Gloves (SMG) certification assurance for quality low-protein latex exam glove.
- Biodegradable and Environment friendly.

| Item | Description | Size | Packaging |
|---------|--------------------------|---------|---------------------------------|
| MG100XS | Medgluv Latex Exam Glove | X-Small | 100 eaches/box 10 boxes/case |
| MG100S | Medgluv Latex Exam Glove | Small | 100 eaches/box 10 boxes/case |
| MG100M | Medgluv Latex Exam Glove | Medium | 100 eaches/box 10 boxes/case |
| MG100L | Medgluv Latex Exam Glove | Large | 100 eaches/box 10 boxes/case |
| MG100XL | Medgluv Latex Exam Glove | X-Large | 100 eaches/box 10 boxes/case |

Product Specifications

| Gauge Thickness | MM | MIL |
|-----------------|-------|------|
| Middle Finger: | .16 | 6.5 |
| Palm: | .13 | 5.2 |
| Cuff: | .10 | 4.0 |
| Average Length | 242mm | 9.5" |

Physical Properties

Before Aging

| | |
|----------------------|--------|
| Tensile Strength: | 29 MPa |
| Ultimate Elongation: | 700% |

After Aging

| | |
|----------------------|--------|
| Tensile Strength: | 25 MPa |
| Ultimate Elongation: | 600% |

Quality Standards

Manufactured in accordance with Quality System ISO 9001.
Meet and exceeds ASTM D3578, D5151 and FDA CFR 177-2600 -
US Food Contact Approved, QSR/FDA 510K Medical Examination Grade.

Note: Specifications are subject to change without notice.



Caution: This product contains natural rubber latex (latex) which may cause allergic reactions in some users. Safe use of this glove by or on latex sensitized individuals has not been established.

For additional information, please visit: www.medgluv.com
Call 1-866-MEDGLUV (1-866-633-4588).

Keep out of sunlight. Store in a cool dry place. Keep away from sources of ozone and ignition.

