

Powder-Free Latex Exam Gloves

NeuGrip Latex Powder-Free Exam Gloves are durable and strong in providing an excellent protection for the aviation and automotive industry. These dependable thicker latex glove provides a secure and excellent grip in any heavy duty environment. They allows easy donning and a comfortable feel and fit. They offer a soft and comfortable feel with outstanding tactile sensitivity and dexterity.

Description:

- NeuGrip's 8 mil glove thickness provides enhanced protection in the finger-tips where protection is critical.
- Unique double chlorination process results in consistency excellent donning, soft, comfortable fit and excellent elasticity
- Standard Malaysian Gloves (SMG) certification assurance for quality low-protein latex exam glove.
- Ideal for tasks in heavy-duty environments..
- Biodegradable and Environment friendly.

| Item | Description | Size | Packaging |
|---------|--------------------------|---------|---------------------------------|
| MG105XS | NeuGrip Latex Exam Glove | X-Small | 100 eaches/box 10 boxes/case |
| MG105S | NeuGrip Latex Exam Glove | Small | 100 eaches/box 10 boxes/case |
| MG105M | NeuGrip Latex Exam Glove | Medium | 100 eaches/box 10 boxes/case |
| MG105L | NeuGrip Latex Exam Glove | Large | 100 eaches/box 10 boxes/case |
| MG105XL | NeuGrip Latex Exam Glove | X-Large | 100 eaches/box 10 boxes/case |

Product Specifications

| Gauge Thickness | MM | MIL |
|-----------------|-------|------|
| Middle Finger: | .20 | 7.9 |
| Palm: | .15 | 6.1 |
| Cuff: | .12 | 4.7 |
| Average Length | 242mm | 9.5" |

Physical Properties

Before Aging

| | |
|----------------------|--------|
| Tensile Strength: | 32 MPa |
| Ultimate Elongation: | 780% |

After Aging

| | |
|----------------------|--------|
| Tensile Strength: | 28 MPa |
| Ultimate Elongation: | 720% |

Quality Standards

Manufactured in accordance with Quality System ISO 9001.
Exceeds and meets current ASTM 3578 standard for critical defects.

Note: Specifications are subject to change without notice.

EXCEPTIONAL FIT,
THICKNESS FOR
DURABLE PROTECTION



Caution: This product contains natural rubber latex which may cause allergic reactions in some users. Safe use of this glove by or on latex sensitized individuals has not been established.

For additional information, please visit: www.medgluv.com
Call 1-866-MEDGLUV (1-866-633-4588).

Keep out of sunlight. Store in a cool dry place. Keep away from sources of ozone and ignition.

