

NeuSkin Vinyl



Powder-Free Vinyl Exam Gloves, Clear

NeuSkin Vinyl exam gloves provide complete hand protection for a variety of applications. This 100% synthetic glove is powder-free and features a smooth surface for more natural tactile sensitivity.

**ENHANCED
FORMULATION
THAT PROVIDES
ELASTICITY AND FIT**

Description:

- NeuSkin powder-free synthetic vinyl glove are a cost effective glove.
- Effective patient care and general usage.
- Low risk situations.
- Smooth fingertips.
- Ambidextrous and powder-free.
- Meets and exceeds ASTM D5250.



Item	Description	Size	Packaging
MG300XS	NeuSkin Clear Vinyl Exam Glove	X-Small	100 eaches/box 10 boxes/case
MG300S	NeuSkin Clear Vinyl Exam Glove	Small	100 eaches/box 10 boxes/case
MG300M	NeuSkin Clear Vinyl Exam Glove	Medium	100 eaches/box 10 boxes/case
MG300L	NeuSkin Clear Vinyl Exam Glove	Large	100 eaches/box 10 boxes/case
MG300XL	NeuSkin Clear Vinyl Exam Glove	X-Large	100 eaches/box 10 boxes/case

Product Specifications

Gauge Thickness	MM	MIL
Middle Finger:	.13	5.2
Palm:	.10	4.0
Cuff:	.08	3.4
Average Length	242mm	9.5"

Physical Properties

Before Aging	
Tensile Strength:	20 MPa
Ultimate Elongation:	450%
After Aging	
Tensile Strength:	19 MPa
Ultimate Elongation:	400%

Quality Standards

Manufactured in accordance with Quality System ISO 9001.
Exceeds current ASTM 5250 Standards for critical defects (AQL 2.5).
AQL for critical defects is 1.5.

Note: Specifications are subject to change without notice.



For additional information, please visit: www.medgluv.com
Call 1-866-MEDGLUV (1-866-633-4588).

Keep out of sunlight. Store in a cool dry place. Keep away from sources of ozone and ignition.



www.medgluv.com

