

# NeuSkin Ultra



## Powder-Free Vinyl Exam Gloves, Stretch

NeuSkin Ultra Stretch Vinyl exam gloves are made from a vinyl formulation that helps provide elasticity and fit. They are soft and flexible, don easily and feel natural. Stretch Vinyl exam gloves are ideal for short-term use where there is minimal stress on the glove and a low-to-no risk of exposure to blood or other potentially infectious materials.

ENHANCED  
FORMULATION  
THAT PROVIDES  
**ELASTICITY  
AND FIT**

### Description:

- NeuSkin Ultra Synthetic vinyl powder-free stretch vinyl exam gloves are an advanced formulation stretch vinyl that conform to the hands for glove that gives outstanding tactile sensitivity with a unique softness and flexibility.
- The stretchy vinyl glove incorporates a unique PVC blend that looks and feels like a latex gloves.
- Tapered at the cuff to minimize cuff roll down. Easy to put on and take off.
- Meets and exceeds ASTM D5250

| Item    | Description                            | Size    | Packaging                       |
|---------|----------------------------------------|---------|---------------------------------|
| MG325L  | NeuSkin Ultra Stretch Vinyl Exam Glove | X-Small | 100 eaches/box<br>10 boxes/case |
| MG325S  | NeuSkin Ultra Stretch Vinyl Exam Glove | Small   | 100 eaches/box<br>10 boxes/case |
| MG325M  | NeuSkin Ultra Stretch Vinyl Exam Glove | Medium  | 100 eaches/box<br>10 boxes/case |
| MG325L  | NeuSkin Ultra Stretch Vinyl Exam Glove | Large   | 100 eaches/box<br>10 boxes/case |
| MG325XL | NeuSkin Ultra Stretch Vinyl Exam Glove | X-Large | 100 eaches/box<br>10 boxes/case |

### Product Specifications

| Gauge Thickness | MM    | MIL  |
|-----------------|-------|------|
| Middle Finger:  | .13   | 5.5  |
| Palm:           | .10   | 4.3  |
| Cuff:           | .08   | 3.4  |
| Average Length  | 242mm | 9.5" |

### Physical Properties

| Before Aging         |        |
|----------------------|--------|
| Tensile Strength:    | 28 MPa |
| Ultimate Elongation: | 450%   |
| After Aging          |        |
| Tensile Strength:    | 22 MPa |
| Ultimate Elongation: | 400%   |

### Quality Standards

Manufactured in accordance with Quality System ISO 9001.  
Exceeds current ASTM 5250 Standards for critical defects (AQL 2.5).  
AQL for critical defects is 1.5.

Note: Specifications are subject to change without notice.



For additional information, please visit: [www.medgluv.com](http://www.medgluv.com)  
Call 1-866-MEDGLUV (1-866-633-4588).

Keep out of sunlight. Store in a cool dry place. Keep away from sources of ozone and ignition.



[www.medgluv.com](http://www.medgluv.com)

