

NitraGrip XP



Extended Cuff PF Nitrile Exam Gloves

NitraGrip XP offer dependable strength you can count on and have a long 12" extended cuff that helps provide extra protection for your wrists and arms. Full textured for enhanced grip. Barrier tested against fentanyl exposure to help protect you when treating overdose patients.

MAXIMUM BARRIER PROTECTION FOR ALL RISKS

Description:

- When a high-risk situation requires the highest standard in protection.
- These elongated synthetic gloves provides extra protection and typically used for tasks involving prolonged exposure to blood, body fluids, surgical instrument cleaning solutions and other chemicals.
- Meets and exceeds NIOSH recommendations for Fentanyl.
- Meets both ASTM D6978-05 and USP 800 Personal Protection Equipment requirement for use of chemotherapy gloves when handling hazardous drugs.



Item	Description	Size	Packaging
MG50051	NitraGrip XP Nitrile Exam Glove	Small	100 eaches/box 10 boxes/case
MG50052	NitraGrip XP Nitrile Exam Glove	Medium	100 eaches/box 10 boxes/case
MG50053	NitraGrip XP Nitrile Exam Glove	Large	100 eaches/box 10 boxes/case
MG50054	NitraGrip XP Nitrile Exam Glove	X-Large	100 eaches/box 10 boxes/case
MG50055	NitraGrip XP Nitrile Exam Glove	2X-Large	90 eaches/box 10 boxes/case
MG50056	NitraGrip XP Nitrile Exam Glove	3X-Large	90 eaches/box 10 boxes/case

Product Specifications

Gauge Thickness	MM	MIL
Middle Finger:	.16	6.3
Palm:	.12	4.7
Cuff:	.10	3.7
Average Length	305mm	12"

Physical Properties

Before Aging

Tensile Strength:	35 MPa
Ultimate Elongation:	600%

After Aging

Tensile Strength:	32 MPa
Ultimate Elongation:	550%

Quality Standards

Meet and exceeds ASTM D6319 and F1671: Standard Specification for Nitrile: Examination Gloves for Medical Application and Standard Test method for Resistance of Materials used in Protective Clothing in Penetration by Blood-Borne Pathogens Using Phi-X174 Bacteriophage as a Test System.

Note: Specifications are subject to change without notice.



Tested with 30 CHEMO DRUGS in accordance with ASTM D6978 and barrier tested against Fentanyl exposure.



This glove meets **NFPA 1999**, Standard on Protective Clothing and Ensembles for Emergency Medical Operations.

For additional information, please visit: www.medgluv.com Call 1-866-MEDGLUV (1-866-633-4588).

Keep out of sunlight. Store in a cool dry place. Keep away from sources of ozone and ignition.



www.medgluv.com

