

## 2-Ply Blue/White pH Nitrile Exam Gloves

These 2-Ply technology outer Blue and inner White elongated latex-free nitrile gloves offers the unique ability to detects tears and rips by showing the high visibility white inner through the breach.

### Description:

- Designed specially for the Emergency Medical Profession for high risk situation that requires the highest standards in double protection and safety in administering critical and patient care.
- Breach alert dual-layer provides extra durability and safety.
- Ideal for EMS, Law Enforcement and Corrections
- Meets and exceeds NIOSH recommendations for Fentanyl and gastric acid resistant.
- Meets both ASTM D6978-05 and USP 800 Personal Protection Equipment requirement for use of chemotherapy gloves when handling hazardous drugs.
- Meets and exceeds the standard emergency glove medical of NFPA 1999 Standards on protective clothing, Emergency Medical Operations, Edition 2018

Item	Description	Size	Packaging
MG50091	NitraGrip Pro, 2-Ply Nitrile Exam Glove	Small	50 eaches/box 10 boxes/case
MG50092	NitraGrip Pro, 2-Ply Nitrile Exam Glove	Medium	50 eaches/box 10 boxes/case
MG50093	NitraGrip Pro, 2-Ply Nitrile Exam Glove	Large	50 eaches/box 10 boxes/case
MG50094	NitraGrip Pro, 2-Ply Nitrile Exam Glove	X-Large	50 eaches/box 10 boxes/case
MG50095	NitraGrip Pro, 2-Ply Nitrile Exam Glove	2X-Large	45 eaches/box 10 boxes/case
MG50096	NitraGrip Pro, 2-Ply Nitrile Exam Glove	3X-Large	45eaches/box 10 boxes/case

### Product Specifications

Gauge Thickness	MM	MIL
Middle Finger:	.21	8.2
Palm:	.17	6.8
Cuff:	.14	5.5
Average Length	305mm	12"

### Physical Properties

Before Aging	
Tensile Strength:	35 MPa
Ultimate Elongation:	600%
After Aging	
Tensile Strength:	32 MPa
Ultimate Elongation:	550%

### Quality Standards

Meet and exceeds ASTM D6319 and F1671: Standard Specification for Nitrile: Examination Gloves for Medical Application and Standard Test method for Resistance of Materials used in Protective Clothing in Penetration by Blood-Borne Pathogens Using Phi-X174 Bacteriophage as a Test System.

Note: Specifications are subject to change without notice.



## EXCELLENT BARRIER PROTECTION FOR ALL RISKS



**Tested with 30 CHEMO DRUGS**  
in accordance with ASTM D6978 and barrier tested against Fentanyl exposure.



This glove meets **NFPA 1999**, Standard on Protective Clothing and Ensembles for Emergency Medical Operations.

For additional information, please visit: [www.medgluv.com](http://www.medgluv.com)  
Call 1-866-MEDGLUV (1-866-633-4588).

Keep out of sunlight. Store in a cool dry place. Keep away from sources of ozone and ignition.

