

NitriSkin White



Powder-Free Nitrile Exam Gloves

NitriSkin White Nitrile are designed mainly for the small facilities physician offices and will offer you the protection of Nitrile and feeling of latex. Featured with an advanced stretch formulation, it will provide you ultimate comfort for extended periods.

Description:

- Textured fingertips enables grip, wet or dry.
- Tear resistant beaded cuff eases donning.
- Standard color coded sizing to identify different sizes.
- Offer excellent sensitivity for touch perception and dependable protection.
- Meets both ASTM D6978-05 and USP 800 Personal Protection Equipment requirement for use of select chemotherapy gloves.

Item	Description	Size	Packaging
MG502XS	NitriSkin White Nitrile Exam Glove	X-Small	100 eaches/box 10 boxes/case
MG502S	NitriSkin White Nitrile Exam Glove	Small	100 eaches/box 10 boxes/case
MG502M	NitriSkin White Nitrile Exam Glove	Medium	100 eaches/box 10 boxes/case
MG502L	NitriSkin White Nitrile Exam Glove	Large	100 eaches/box 10 boxes/case
MG502XL	NitriSkin White Nitrile Exam Glove	X-Large	100 eaches/box 10 boxes/case

Product Specifications

Gauge Thickness	MM	MIL
Middle Finger:	.09	3.5
Palm:	.08	3.1
Cuff:	.06	2.4
Average Length	242mm	9.5"

Physical Properties

Before Aging	
Tensile Strength:	38 MPa
Ultimate Elongation:	580%
After Aging	
Tensile Strength:	35 MPa
Ultimate Elongation:	550%


Quality Standards

Exceed current ASTM D6319: Standard Specification for Nitrile: Examination Gloves for Medical Application.
No evidence of dermal irritation or contact sensitization as measured by Primary Skin Irritation.

Note: Specifications are subject to change without notice.

EXCELLENT
BARRIER PROTECTION
FOR ALL RISKS



 **Tested with 30 CHEMO DRUGS**
in accordance with ASTM D6978 and barrier tested against Fentanyl exposure.

For additional information, please visit: www.medgluv.com
Call 1-866-MEDGLUV (1-866-633-4588).

Keep out of sunlight. Store in a cool dry place. Keep away from sources of ozone and ignition.

