

# OysterSkin



## Powder-Free Nitrile Exam Gloves

OysterSkin Grey Nitrile Exam Gloves provide tactile sensitivity and protects hands from potentially infectious materials. New nitrile compounding technology allows for a thinner glove with incredible elasticity. They offer a soft and comfortable feel with outstanding tactile sensitivity. Packaged in 250 gloves per box. Devotes less space to glove storage, less time handling cartons

### Description:

- Finger-tip textured enable grip even when wet.
- Virtually no pressure on the hands and reduces fatigue after extended wear.
- Offer excellent sensitivity for touch perception and dependable protection.
- 60% saving in the storage space, fewer boxes to discard and less waste.
- Meets both ASTM D6978-05 and USP 800 Personal Protection Equipment requirement for use of chemotherapy gloves when handling hazardous drugs.
- Meets and exceeds NIOSH recommendations for Fentanyl.

Item	Description	Size	Packaging
MG5301	OysterSkin Grey Nitrile Exam Glove	X-Small	250 eaches/box 10 boxes/case
MG5302	OysterSkin Grey Nitrile Exam Glove	Small	250 eaches/box 10 boxes/case
MG5303	OysterSkin Grey Nitrile Exam Glove	Medium	250 eaches/box 10 boxes/case
MG5304	OysterSkin Grey Nitrile Exam Glove	Large	250 eaches/box 10 boxes/case
MG5305	OysterSkin Grey Nitrile Exam Glove	X-Large	230 eaches/box 10 boxes/case

### Product Specifications

Gauge Thickness	MM	MIL
Middle Finger:	.09	3.5
Palm:	.08	3.1
Cuff:	.06	2.4
Average Length	242mm	9.5"

### Physical Properties

#### Before Aging

Tensile Strength:	38 MPa
Ultimate Elongation:	580%

#### After Aging

Tensile Strength:	35 MPa
Ultimate Elongation:	550%

### Quality Standards

Exceed current ASTM D6319: Standard Specification for Nitrile: Examination Gloves for Medical Application.

No evidence of dermal irritation or contact sensitization as measured by Primary Skin Irritation.

Note: Specifications are subject to change without notice.

## EXCELLENT BARRIER PROTECTION FOR ALL RISKS



**Tested with 30 CHEMO DRUGS**  
in accordance with ASTM D6978 and barrier tested against Fentanyl exposure.

For additional information, please visit: [www.medgluv.com](http://www.medgluv.com)  
Call 1-866-MEDGLUV (1-866-633-4588).

Keep out of sunlight. Store in a cool dry place. Keep away from sources of ozone and ignition.



[www.medgluv.com](http://www.medgluv.com)

