

Powder-Free Nitrile Exam Gloves

NitraCare Nitrile Exam Gloves are our best-selling exam gloves in our full exam glove portfolio. They provide excellent barrier protection from bacteria, viruses and chemicals during procedures where risk of fluid exposure is moderate to high.

Description:

- Textured fingertips enables grip, wet or dry.
- Tear resistant beaded cuff eases donning.
- Standard color coded sizing to identify different sizes.
- Offer excellent sensitivity for touch perception and dependable protection.
- Meets both ASTM D6978-05 and USP 800 Personal Protection Equipment requirement for use of select chemotherapy gloves.

| Item | Description | Size | Packaging |
|---------|------------------------------|----------|---------------------------------|
| MNE5051 | NitraCare Nitrile Exam Glove | X-Small | 100 eaches/box 10 boxes/case |
| MNE5052 | NitraCare Nitrile Exam Glove | Small | 100 eaches/box 10 boxes/case |
| MNE5053 | NitraCare Nitrile Exam Glove | Medium | 100 eaches/box 10 boxes/case |
| MNE5054 | NitraCare Nitrile Exam Glove | Large | 100 eaches/box 10 boxes/case |
| MNE5055 | NitraCare Nitrile Exam Glove | X-Large | 100 eaches/box 10 boxes/case |
| MNE5056 | NitraCare Nitrile Exam Glove | 2X-Large | 90 eaches/box 10 boxes/case |

Product Specifications

| Gauge Thickness | MM | MIL |
|-----------------|-------|------|
| Middle Finger: | .09 | 3.5 |
| Palm: | .08 | 3.1 |
| Cuff: | .06 | 2.4 |
| Average Length | 242mm | 9.5" |

Physical Properties

Before Aging

| | |
|----------------------|--------|
| Tensile Strength: | 35 MPa |
| Ultimate Elongation: | 580% |

After Aging

| | |
|----------------------|--------|
| Tensile Strength: | 38 MPa |
| Ultimate Elongation: | 550% |

Quality Standards

Exceed current ASTM D6319: Standard Specification for Nitrile: Examination Gloves for Medical Application.

No evidence of dermal irritation or contact sensitization as measured by Primary Skin Irritation.

Note: Specifications are subject to change without notice.



EXCELLENT BARRIER PROTECTION FOR ALL RISKS



Tested with 30 CHEMO DRUGS
in accordance with ASTM D6978 and barrier tested against Fentanyl exposure.

For additional information, please visit: www.medgluv.com
Call 1-866-MEDGLUV (1-866-633-4588).

Keep out of sunlight. Store in a cool dry place. Keep away from sources of ozone and ignition.

