

NitriSkin Sterile Singles and Pairs

Powder-Free Sterile Examination Gloves

NitriSkin Sterile Singles and Pairs are now the ultimate choice use of medical professionals in high demand environments like Obstetricians and Gynecologists, intensive care unit ICC/CCU, traumas, burn centers and wound care that requires a high quality latex-free sterile glove that delivers maximum barrier protection.

Characteristics and Overview:

- Ideal for sterile clinical procedures when an exceptionally resistant non-latex glove is required.
- Textured gloves enhance grip in wet and dry environments.
- Exceptionally stretchy soft, low modulus nitrile synthetic compound is very comfortable for extended periods.
- Compact packaging makes these sterile gloves convenient enough to place in a pocket and easy to store.
- Meets ASTM D6978-05 for use of chemotherapy gloves when handling hazardous drugs.



Item No.	Description	Size	Packaging
MGSE502S	NitriSkin Sterile Singles Nitrile Exam Glove	Small	100 eaches/box 4 boxes/case
MGSE502M	NitriSkin Sterile Singles Nitrile Exam Glove	Medium	100 eaches/box 4 boxes/case
MGSE502L	NitriSkin Sterile Singles Nitrile Exam Glove	Large	100 eaches/box 4 boxes/case
MGSE502XL	NitriSkin Sterile Singles Nitrile Exam Glove	X-Large	100 eaches/box 4 boxes/case
MGSE552S	NitriSkin Sterile Pairs Nitrile Exam Glove	Small	50 pairs/box 4 boxes/case
MGSE552M	NitriSkin Sterile Pairs Nitrile Exam Glove	Medium	50 pairs/box 4 boxes/case
MGSE552L	NitriSkin Sterile Pairs Nitrile Exam Glove	Large	50 pairs/box 4 boxes/case
MGSE552XL	NitriSkin Sterile Pairs Nitrile Exam Glove	X-Large	50 pairs/box 4 boxes/case



Glove Specifications:

Brand	NitriSkin
Material	Nitrile
Color	Sky Blue
Sterility	Sterile
Cuff	Beaded Cuff
Sizes	SM, MD, LG and XL
Accelerators	No
Chemo Tested	Yes
Glove Interior	Powder-Free
Glove Exterior	Full Textured
Length (inches)	9.5 in
Palm Thickness (Avg)	4.3 mil
Finger Thickness (Avg)	5.1 mil
Tensile Strength (Mpa)	35 / 33 Mpa
Ultimate Elongation (%)	600 / 560
Pin Hole (AQL)	1.5

Quality Standards:

Manufactured in accordance with
Quality System ISO 9001.

Exceeds and meets current ASTM Standards:
- D6319, D5151, D1671 & D6978.

*Summary of Chemotherapy Testing Reports - ASTM D6978-05

Drug Name	Average Breakthrough Detection Time
Carboplatin	No breakthrough up to 240 min
Cisplatin	No breakthrough up to 240 min
Cyclophosphamide	No breakthrough up to 240 min
Dacarbazine	No breakthrough up to 240 min
Doxorubicin	No breakthrough up to 240 min
Etoposide	No breakthrough up to 240 min
Fluorouracil	No breakthrough up to 240 min
Ifosfamide	No breakthrough up to 240 min
Methotrexate	No breakthrough up to 240 min
Mitomycin C	No breakthrough up to 240 min
Mitoxantrone	No breakthrough up to 240 min
Paclitaxel	No breakthrough up to 240 min
Vincristine Sulfate	No breakthrough up to 240 min
Carmustine	Not for use with Carmustine
ThioTepa	Not for use with ThioTepa