

NitraGrip Sterile Pairs

12" Powder-Free Sterile Examination Gloves

NitraGrip Sterile glove is made from a synthetic, nitrile formulation that provides high elasticity with superior barrier protection and puncture resistance. It is an ideal solution for a range of sterile examination (non-surgical) procedures.

Characteristics and Overview:

- Perfect for USP 797 and 800 standards compliance for pharmacy compounding and sterile processing.
- Sterile procedures outside of the O.R.
- Meets and exceeds NIOSH recommendations for Fentanyl.
- Meets both ASTM D6978-05 and USP 800 Standards for use of chemotherapy gloves when handling hazardous drugs*



Item No.	Description	Size	Packaging
MGSE5121	NitraGrip Sterile Pairs Nitrile Exam Glove	Small	50 pairs/box 4 boxes/case
MGSE5122	NitraGrip Sterile Pairs Nitrile Exam Glove	Medium	50 pairs/box 4 boxes/case
MGSE5123	NitraGrip Sterile Pairs Nitrile Exam Glove	Large	50 pairs/box 4 boxes/case
MGSE5124	NitraGrip Sterile Pairs Nitrile Exam Glove	X-Large	50 pairs/box 4 boxes/case



Glove Specifications:

Brand	NitraGrip
Material	Nitrile
Color	White
Sterility	Non-Sterile
Cuff	Beaded Cuff
Sizes	SM, MD, LG and XL
Accelerators	Yes
Chemo Tested	Yes
Glove Interior	Powder-Free
Glove Exterior	Full Textured
Length (inches)	12 in
Palm Thickness (Avg)	4.5 mil
Finger Thickness (Avg)	6.3 mil
Tensile Strength (Mpa)	35 / 33 Mpa
Ultimate Elongation (%)	600 / 560
Pin Hole (AQL)	1.5

*Summary of Chemotherapy Testing Reports - ASTM D6978-05

Drug Name	Average Breakthrough Detection Time
Carboplatin	No breakthrough up to 240 min
Cisplatin	No breakthrough up to 240 min
Cyclophosphamide	No breakthrough up to 240 min
Dacarbazine	No breakthrough up to 240 min
Doxorubicin	No breakthrough up to 240 min
Etoposide	No breakthrough up to 240 min
Fluorouracil	No breakthrough up to 240 min
Ifosfamide	No breakthrough up to 240 min
Methotrexate	No breakthrough up to 240 min
Mitomycin C	No breakthrough up to 240 min
Mitoxantrone	No breakthrough up to 240 min
Paclitaxel	No breakthrough up to 240 min
Vincristine Sulfate	No breakthrough up to 240 min
Carmustine	Not for use with Carmustine
ThioTepa	Not for use with ThioTepa

Quality Standards:

Manufactured in accordance with Quality System ISO 9001.
 Exceeds and meets current ASTM Standards:
 - D6319, D5151, D1671 & D6978.