

## NeuMask

### Ear-loop Procedure Mask

Medgluv Procedure Earloop Face Masks are a level 3 high barrier mask ideal for oral exams of high-risk patients and procedures involving high levels of spray, spatter, moisture and airborne particles. Perfect for ultrasonic scalers, high-speed aerosols, air polishers and oral surgeries.

Soft and strong, ample-length earloops so there is less pulling, and adjustable nosepiece for comfortable fit. Especially useful for laser procedures that generate harmful airborne particles. Comfortable and highly breathable in any settings.

Light, comfortable masks for extended wear when necessary. Latex free, ear loops are made with no irritating rubber or plastic touches the skin.

Three-ply, glass-free filter media is surrounded by soft, non-irritating inner and outer facing.

### Description

- Procedure masks with elastic ear-loops provide complete protection
- Polypropylene material offer a soft, fluid protection barrier.
- A full-width nosepiece guarantees proper fit.
- High Filtration Melt-blown Filter. Easy breathing; Delta-P < 6.0
- Bacterial filtration Efficiency (BFE) > 98% @ 3.0µm.
- Meets and exceeds ASTM F2100-19 Standards

Item No: MGM101 (Level 3)  
 Color: Blue  
 Packaging: 50 masks/box, 10 boxes/case



Item No: MGM102 (Level 2)  
 Color: Orange  
 Packaging: 50 masks/box, 10 boxes/case



Item No: MGM103 (Level 2)  
 Color: Pink  
 Packaging: 50 masks/box, 10 boxes/case



Item No: MGM104 (Level 3)  
 Color: Black  
 Packaging: 50 masks/box, 10 boxes/case



Item No: MGMC100 (Level 1)  
 Color: Childs Characters  
 Packaging: 50 masks/box, 10 boxes/case

