

NeuMask Procedure Ear-loop Facemask

Medgluv NeuMask Procedure Ear-loop face masks are comfortable and highly breathable. This product is appropriate for all general patient care situations.

Characteristics and Overview:

- Basic Procedure masks with elastic ear-loops provide complete protection while remaining light and comfortable.
- Polypropylene outer facing on some procedure masks offer a soft, fluid protection barrier while the inner facing assures comfort and breathability.
- A full-width nosepiece guarantees proper fit.
- High Filtration Meltblown Filter
- Bacterial filtration Efficiency (BFE) > 98% @ 3.0µm
- Easy breathing; Delta-P < 6.0
- Meets and exceeds ASTM F2100-19 Level 2 Standards.



Mask Specifications:

Brand	NeuMask
Material	Polypropylene
Color	Blue
Packaging	50 masks/box, 10 boxes/case



SureFit Glove Dispenser

Characteristics and Overview:

- White-coated wire glove box holders
- Horizontal mount, can fit one or two boxes.

Dispenser Specifications:

Brand	SureFit
Material	Wire
Color	Blue
Depth	4 inches
Width	10.5 inches
Packaging	10 dispensers/box

