

Medgluv Shoe Covers

With Non-Skid Bottom

Medgluv disposable shoe covers are a great way to keep dirt and other contaminants out of unwanted areas in hospitals, institutions, labs, clean room or other manufacturing facilities. Prevent slips and contamination, applicable in a wide variety of settings and situations.

Characteristics and Overview:

- Increased tear strength with sewn seams.
- Skid-resistant tacky bottom strip
- Full elastic added to fit over shoes for easy donning.
- Sewn Seams for Superior Strength.
- Bacterial filtration Efficiency (BFE) > 98% @ 3.0µm
- Non-Conductive Material
- Latex-Free



Mask Specifications:


Brand:	Medgluv Shoe Covers
Material:	Spunbond non-woven material
Color:	Blue
Item Number:	MGSC102 (Size XL)
Packaging:	100covers/bag, 10 bags/case

Medgluv Bouffant Caps

Medgluv Bouffant Caps are made with lightweight spunbond fabric and are designed to cover hair in no fluid contamination or light fluid contamination environments. Regular blue spunbond bouffants, are suitable for non-fluid contamination environments, such as patient rooms, waiting areas and exam rooms. They are roomy and lightweight to wear, keeping hair securely tucked away.

Characteristics and Overview:

- One layer spunbond fabric
- Latex-Free elastic
- Blue Color
- No potential fluid contact.



Caps Specifications:

Brand	Medgluv Bouffant Caps
Material	Spunbond non-woven material
Color	Blue
Item Number:	MGBC21 - 21" MGMC24 - 24"
Packaging	100caps/bag, 10 bags/case



Single use only
Do not reuse.



Latex-Free



Non-Sterile