

Powder-Free Nitrile Exam Gloves



Diamond Textured

SonicBurst are perfect for EMS, First Responders and Law Enforcements. This glove offers dependable strength and diamond textured grip you can count on providing extra protection. Barrier tested against fentanyl exposure to help protect you when treating overdose patients.

Description:

- When a high-risk situation requires the highest standard in protection.
- Offer excellent sensitivity for touch perception and dependable protection
- Meet standards of the National Fire Protection Association (NFPA)
- Meets and exceeds NIOSH recommendations for Fentanyl.
- Meets both ASTM D6978-05 and USP 800 Personal Protection Equipment requirement for use of chemotherapy gloves when handling hazardous drugs.

Item	Description	Size	Packaging
MG5070	SonicBurst Nitrile Exam Glove	X-Small	100 eaches/box 10 boxes/case
MG5071	SonicBurst Nitrile Exam Glove	Small	100 eaches/box 10 boxes/case
MG5072	SonicBurst Nitrile Exam Glove	Medium	100 eaches/box 10 boxes/case
MG5073	SonicBurst Nitrile Exam Glove	Large	100 eaches/box 10 boxes/case
MG5074	SonicBurst Nitrile Exam Glove	X-Large	100 eaches/box 10 boxes/case
MG5075	SonicBurst Nitrile Exam Glove	2X-Large	90 eaches/box 10 boxes/case
MG5076	SonicBurst Nitrile Exam Glove	3X-Large	90 eaches/box 10 boxes/case



Product Specifications

Gauge Thickness	MM	MIL
Middle Finger:	.21	6.5
Palm:	.17	5.5
Cuff:	.14	4.9
Average Length	230mm	9.5"

Physical Properties

Before Aging

Tensile Strength:	38 MPa
Ultimate Elongation:	650%

After Aging

Tensile Strength:	35 MPa
Ultimate Elongation:	600%

Quality Standards

Meet and exceeds ASTM D6319 and F1671: Standard Specification for Nitrile Examination Gloves for Medical Application and Standard Test method for Resistance of Materials used in Protective Clothing in Penetration by Blood-Borne Pathogens Using Phi-X174 Bacteriophage as a Test System.

Note: Specifications are subject to change without notice.



Tested with 30 CHEMO DRUGS

in accordance with ASTM D6978 and barrier tested against Fentanyl exposure.



This glove meets **NFPA 1999**, Standard on Protective Clothing and Ensembles for Emergency Medical Operations.

For additional information, please visit: www.medgluv.com
Call 1-866-MEDGLUV (1-866-633-4588).

Keep out of sunlight. Store in a cool dry place. Keep away from sources of ozone and ignition.

