

Medgluv Child Face Mask

Ear-loop Procedure Mask, Level 1

KIDS MASK THAT MEETS FDA GUIDELINES FOR PEDIATRIC USE

The CDC and Prevention protocols for respiratory etiquette recommends that kids with undiagnosed respiratory symptoms wear face masks in waiting rooms. Any child to hospital admission in waiting areas should have a child-size mask available for kids. Medgluv pediatric mask meets FDA Guidelines for use with children. Medgluv new latex-free pediatric mask helps protect patients from colds, flu and other respiratory threats.



KID'S MASKS WITH HAPPY FUN CHARACTERS

Description

- Procedure masks with elastic ear-loops provide complete protection
- Soft lightweight materials for comfort wear
- Nosepiece for custom fit on small faces
- Meets and exceeds ASTM F2100-19 Level 1 Standards.

Item No: MGMC100 (Level 1)

Color: Kids Characters

Packaging: 50 masks/box
10 boxes/case

** Bacterial filtration Efficiency (BFE) > 96% @ 3.0µm.

