

NeuPro 16



Elongated Nitrile Exam Gloves

NeuPro 16 gloves are specially designed for sterilization, decontamination, disinfection, cleaning, and handling chemicals in hazardous environments. This 16-inch elongated specialty glove provides exceptional extended coverage and protection, particularly in situations with potential exposure to sharps and harmful chemicals.


Tested against a wide range of chemotherapy drugs, NeuPro 16 gloves offer reliable defense and peace of mind when handling instruments and decontamination agents.

THE NEW STANDARDS FOR PROTECTION IN THE STERILE PROCESSING DEPARTMENT.

Description:

- Palm textured and finger textured to enhance grip in wet and dry situations.
- Longer 16" extended length
- Tested for 48 chemotherapy drugs.
- Tested for puncture resistance and color-coded sizing on box.
- Meets and exceeds NIOSH recommendations for Fentanyl and gastric acid resistant.
- Meets both ASTM D6978-05 and USP 800 Personal Protection Equipment requirement for use of chemotherapy gloves when handling hazardous drugs.

Item	Description	Size	Packaging
MG50171	NeuPro 16 Nitrile Exam Glove	Small	50 eaches/box 10 boxes/case
MG50172	NeuPro 16 Nitrile Exam Glove	Medium	50 eaches/box 10 boxes/case
MG50173	NeuPro 16 Nitrile Exam Glove	Large	50 eaches/box 10 boxes/case
MG50174	NeuPro 16 Nitrile Exam Glove	X-Large	50 eaches/box 10 boxes/case



Product Specifications

Gauge Thickness	MM	MIL
Middle Finger:	0.28	11.0
Palm:	0.251	9.9
Cuff:	0.16	6.0
Average Length	410 min	16"

Physical Properties

Before Aging

Tensile Strength:	20 MPa
Ultimate Elongation:	650%

After Aging

Tensile Strength:	16 MPa
Ultimate Elongation:	550%

Quality Standards

Manufactured in compliance with multiple international standards, including ASTM D6319, ASTM F1671, ASTM D5151 and ASTM D6978

Note: Specifications are subject to change without notice.

view this product
on our website:



Tested with 48 CHEMO DRUGS
in accordance with ASTM D6978 and barrier
tested against Fentanyl exposure.

For additional information, please visit: www.medgluv.com
Call 1-866-MEDGLUV (1-866-633-4588).

Keep out of sunlight. Store in a cool dry place. Keep away from sources of ozone and ignition.



www.medgluv.com

