

NeuMask

Ear-loop Procedure Mask, Level 3

Medgluv Procedure Earloop Face Masks are a level 3 high barrier mask ideal for oral exams of high-risk patients and procedures involving high levels of spray, spatter, moisture and airborne particles. Perfect for ultrasonic scalers, high-speed aerosols, air polishers and oral surgeries.

Soft and strong, ample-length earloops so there is less pulling, and adjustable nosepiece for comfortable fit. Especially useful for laser procedures that generate harmful airborne particles. Comfortable and highly breathable in any settings.

Light, comfortable masks for extended wear when necessary. Latex free, ear loops are made with no irritating rubber or plastic touches the skin.

Three-ply, glass-free filter media is surrounded by soft, non-irritating inner and outer facing



Made in USA



Description

- Procedure masks with elastic ear-loops provide complete protection
- Plastic laminated nose wire.
- A full-width nosepiece guarantees proper fit.
- High Filtration Melt-blown Filter. Easy breathing; Delta-P < 6.0
- Bacterial filtration Efficiency (BFE) > 98% @ 3.0µm.
- Meets and exceeds ASTM F2100-23 Level 3

Item No: MGM101 (Level 3)
Color: Blue
Packaging: 50 masks/box, 10 boxes/case



Performance Metrics	Level 1 Mask	Level 2 Mask	Level 3 Mask
PFE at 0.1 micron	≥ 95%	≥ 98%	≥ 98%
BFE at 3.0 microns	≥ 95%	≥ 98%	≥ 98%
Delta P, mm H₂O/cm²	< 5.0	< 6.0	< 6.0
Fluid Resistance, mmHg	≥ 80 mm Hg	≥ 120 mm Hg	≥ 160 mm Hg
Flame Speed	Class 1	Class 1	Class 1

