

NeuArmor



Powder-Free Nitrile Exam Gloves

NeuArmor black nitrile glove is a high performing, latex-free exam glove for use in every application in every environment. Exceptionally stretchy soft, low modulus nitrile synthetic compound excellent for extended periods. Slightly tacky exterior finish provides excellent gripping properties.

Description:

- NeuArmor black exam gloves are extra thicker gloves that provide the advanced protection and optimal grip needed for demanding, high-stress situations.
- Powder-Free nitrile gloves offer dependable strength and are raised textured for an enhanced grip in demanding environments.
- Meets and exceeds NIOSH recommendations for Fentanyl.
- Meets both ASTM D6978-05 and USP 800 Personal Protection Equipment requirement for use of chemotherapy gloves when handling hazardous drugs.

Item	Description	Size	Packaging
MG5901	NeuArmor Black Nitrile Exam Glove	Small	100 eaches/box 10 boxes/case
MG5902	NeuArmor Black Nitrile Exam Glove	Medium	100 eaches/box 10 boxes/case
MG5903	NeuArmor Black Nitrile Exam Glove	Large	100 eaches/box 10 boxes/case
MG5904	NeuArmor Black Nitrile Exam Glove	X-Large	100 eaches/box 10 boxes/case
MG5905	NeuArmor Black Nitrile Exam Glove	XX-Large	90 eaches/box 10 boxes/case

Product Specifications

Gauge Thickness	MM	MIL
Middle Finger:	.14	8.2
Palm:	.12	7.4
Cuff:	.09	6.5
Average Length	242mm	9.5"

Physical Properties


Before Aging	
Tensile Strength:	25 MPa
Ultimate Elongation:	550%
After Aging	
Tensile Strength:	25 MPa
Ultimate Elongation:	500%

Quality Standards

THIS GLOVE IS FOR SINGLE USE ONLY. THIS GLOVE MEETS THE SINGLE-USE EMERGENCY MEDICAL EXAMINATION GLOVE REQUIREMENTS OF NFPA 1999, STANDARD ON PROTECTIVE CLOTHING AND ENSEMBLES FOR EMERGENCY MEDICAL OPERATIONS, 2018 EDITION.

Note: Specifications are subject to change without notice.

MAXIMUM BARRIER PROTECTION FOR ALL RISKS



view this product on our website:



National Fire Protection Agency.



Tested with 48 CHEMO DRUGS
in accordance with ASTM D6978 and barrier tested against Fentanyl exposure.

For additional information, please visit: www.medgluv.com
Call 1-866-MEDGLUV (1-866-633-4588).

Keep out of sunlight. Store in a cool dry place. Keep away from sources of ozone and ignition.



www.medgluv.com

